

Lateral Flow Assay Uncut Sheets

Product Catalogue

Rapid Test

Series	Cat. No.	Uncut Sheet	Specimen	Detection Time
Infectious Diseases	SS801-1	Anti-HIV (1/2) Rapid Test Uncut Sheet	W/S/P	15 min
	SS801-4	HIV Rapid Test Uncut Sheet (4th Gen)	W/S/P	15 min
	SS802-1	HBsAg Rapid Test Uncut Sheet	W/S/P	15 min
	SS803-1	Anti-HCV Rapid Test Uncut Sheet	W/S/P	15 min
	SS804-1	Anti-TP Rapid Test Uncut Sheet	W/S/P	15 min
Respiratory Diseases	SS811-1	SARS-CoV-2 Rapid Test Uncut Sheet	Nasal/Throat swab/Saliva	10-15 min
	SS811-2	SARS-CoV-2 Rapid Test Uncut Sheet	Nasal/Throat swab/Saliva	10-15 min
	SS832-1	Influenza A&B SARS-CoV-2 Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
	SS832-4	RSV/ADV Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
	SS808-1	MP Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
	SS808-2	Anti-MP Rapid Test Uncut Sheet	W/S/P	15 min
	SS809-1	Flu A+B Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
	SS812-1	RSV Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
	SS831-1	Strep A Rapid Test Uncut Sheet	Throat swab	10 min
	SS818-1	ADV Rapid Test Uncut Sheet	Nasal/Throat swab	10-15 min
Tropical Diseases	SS810-1	Malaria P.f Rapid Test Uncut Sheet	Whole blood	15 min
	SS810-2	Malaria (P.f/P.v) Rapid Test Uncut Sheet	Whole blood	15 min
	SS810-3	Malaria (P.f/Pan) Rapid Test uncut sheet	Whole blood	15 min
	SS810-4	Malaria p.f p.LDH Rapid Test uncut sheet	Whole blood	15 min
	SS810-5	Malaria p.f HRPII/p.f pLDH Rapid Test uncut sheet	Whole blood	15 min
	SS806-2	Anti-Dengue IgG/IgM Rapid Test Uncut Sheet	W/S/P	15 min
	SS806-1	Dengue NS1 Rapid Test Uncut Sheet	W/S/P	15 min
Gastrointestinal Diseases	SS805-2	Anti-HP Ab Rapid Test Uncut Sheet	W/S/P	15 min
	SS231-1	FOB Rapid Test Uncut Sheet	Stool	15 min
Cardiac Markers	SS101-1	cTnI Rapid Test Uncut Sheet	W/S/P	15 min

FIA

Series	Cat. No.	Uncut Sheet	Specimen	Linear Range	Detection Time
Respiratory Diseases	SS811-4	SARS-CoV-2 FIA Uncut Sheet	Nasal/Throat swab/Saliva	Qualitative	10-15 min
	SS809-8	Influenza A&B SARS-CoV-2 FIA Uncut Sheet	Nasal/Throat swab	Qualitative	10-15 min
	SS812-3	Flu A+B FIA Uncut Sheet	Nasal/Throat swab/Saliva	Qualitative	10-15 min
	SS818-3	RSV FIA Uncut Sheet	Nasal/Throat swab	Qualitative	10-15 min
	SS831-3	Strep A FIA Uncut Sheet	Throat swab	Qualitative	10 min
	SS811-4	ADV FIA Uncut Sheet	Nasal/Throat swab	Qualitative	10-15 min

Series	Cat. No.	Uncut Sheet	Specimen	Linear Range	Detection Time
Cardiac Markers	SS101-2	cTnI FIA Uncut Sheet	W/S/P	0.1-40 ng/mL	15 min
	SS101-3	hs-cTnI FIA Uncut Sheet	W/S/P	10-10000 pg/mL	15 min
	SS121-1	Myo FIA Uncut Sheet	W/S/P	5-500 ng/mL	15 min
	SS120-1	CK-MB FIA Uncut Sheet	W/S/P	1-100 ng/mL	15 min
	SS230-1	Hcy FIA Uncut Sheet	W/S/P	3.5-48.0 μ mol/L	10 min
	SS102-1	NT-proBNP FIA Uncut Sheet	W/S/P	50-35000 pg/mL	15 min
	SS105-1	D-Dimer FIA Uncut Sheet	W/P	0.1-10 mg/L	5 min
	SS257-1	cTnI/CK-MB/Myo FIA Uncut Sheet	W/S/P	cTnI: 0.1-40 ng/mL Myo: 5-500 ng/mL CK-MB: 1-100 ng/mL	15 min
	SS257-2	cTnI/D-Dimer/NT-proBNP FIA Uncut Sheet	W/P	cTnI: 0.1-40 ng/mL D-Dimer: 0.1-10 mg/L NT-proBNP: 50-35000 pg/mL	15 min
Inflammatory Markers	SS103-1	PCT FIA Uncut Sheet	W/S/P	0.1-100 ng/mL	15 min
	SS122-1	IL-6 FIA Uncut Sheet	W/S/P	5-5000 pg/mL	15 min
	SS106-1	CRP FIA Uncut Sheet	W/S/P	0.5-200 mg/L	3 min
	SS141-1	SAA FIA Uncut Sheet	W/S/P	1-300 mg/L	5 min
Diabetes	SS149-1	HbA1c FIA Uncut Sheet	Whole blood	3-14%	5 min
	SS127-1	C-Peptide FIA Uncut Sheet	W/S/P	0.1-30 ng/mL	15 min
	SS127-2	Insulin FIA Uncut Sheet	W/S/P	2-300 μ U/mL	10 min
Bone Metabolism	SS129-2	25-OH-VD FIA Uncut Sheet (Sandwich Method)	W/S/P	5-100 ng/mL	15 min
Anemia	SS170-1	FERR FIA Uncut Sheet	W/S/P	5-500 ng/mL	10 min
	SS198-1	VB12 FIA Uncut Sheet	S/P	100-2000 pg/mL	10 min

FIA

Series	Cat. No.	Uncut Sheet	Specimen	Linear Range	Detection Time
Preeclampsia	SS131-1	PLGF FIA Uncut Sheet	W/S/P	10-10000 pg/mL	15 min
	SS132-1	sFlt-1 FIA Uncut Sheet	W/S/P	100-50000 pg/mL	15 min
Tumor Markers	SS175-1	CEA FIA Uncut Sheet	W/S/P	1-500 ng/mL	15 min
	SS114-2	fPSA FIA Uncut Sheet	W/S/P	0.2-30 ng/mL	15 min
	SS114-1	PSA FIA Uncut Sheet	W/S/P	2-80 ng/mL	15 min
	SS174-1	AFP FIA Uncut Sheet	W/S/P	5-400 ng/mL	15 min
	SS107-1	CA125 FIA Uncut Sheet	W/S/P	5-3000 U/mL	15 min
Sex Hormones	SS115-1	AMH FIA Uncut Sheet	W/S/P	0.1-20 ng/mL	15 min
	SS157-1	LH FIA Uncut Sheet	W/S/P	1-100 mIU/mL	15 min
	SS156-1	FSH FIA Uncut Sheet	W/S/P	1-100 mIU/mL	15 min
	SS153-1	PRL FIA Uncut Sheet	W/S/P	1-200 ng/mL	15 min
	SS125-1	TESTO FIA Uncut Sheet	W/S/P	0.2-15 ng/mL	10 min
	SS124-1	PROG FIA Uncut Sheet	W/S/P	0.5-60 ng/mL	10 min
	SS148-1	β-HCG FIA Uncut Sheet	W/S/P	5-200000 mIU/mL	15 min
Thyroid	SS147-1	TSH FIA Uncut Sheet	W/S/P	0.1-100 mIU/L	15 min
	SS123-1	T3 FIA Uncut Sheet	S/P	0.61-9.22 nmol/L	10 min
	SS123-3	T4 FIA Uncut Sheet	S/P	12.87-310 nmol/L	15 min
	SS123-2	FT3 FIA Uncut Sheet	S/P	0.60-45 pmol/L	10 min
	SS123-4	FT4 FIA Uncut Sheet	S/P	1-100 pmol/L	15 min
	SS123-5	T3 FIA Uncut Sheet (Sandwich Method)	S/P	0.61-9.22 nmol/L	15 min
	SS123-6	T4 FIA Uncut Sheet (Sandwich Method)	S/P	12.87-310 nmol/L	15 min
	SS123-7	FT3 FIA Uncut Sheet (Sandwich Method)	S/P	0.60-45 pmol/L	15 min
	SS123-8	FT4 FIA Uncut Sheet (Sandwich Method)	S/P	1-100 pmol/L	15 min
Others	SS128-1	PGI/PGII FIA Uncut Sheet	W/S/P	PGI: 5-160 ng/mL PGII: 2-60 ng/mL	15 min
	SS213-1	Total IgE FIA Uncut Sheet	W/S/P	5-5000 ng/mL	15 min

Drug Abuse

Cat. No.	Uncut Sheet	Specimen	Specification	Cut Off
SS189-1-U	MOP Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS189-2-U	MOP Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS189-1-S	MOP Rapid Test Uncut Sheet	Saliva	60*300 mm	40 ng/mL
SS177-1-U	MET Rapid Test Uncut Sheet	Urine	60*300 mm	1000 ng/mL
SS177-2-U	MET Rapid Test Uncut Sheet	Urine	78*300 mm	1000 ng/mL
SS177-1-S	MET Rapid Test Uncut Sheet	Saliva	60*300 mm	25 ng/mL
SS243-1-U	KET Rapid Test Uncut Sheet	Urine	60*300 mm	1000 ng/mL
SS243-2-U	KET Rapid Test Uncut Sheet	Urine	78*300 mm	1000 ng/mL
SS243-1-S	KET Rapid Test Uncut Sheet	Saliva	60*300 mm	100 ng/mL
SS181-1-U	THC Rapid Test Uncut Sheet	Urine	60*300 mm	20 ng/mL
SS181-2-U	THC Rapid Test Uncut Sheet	Urine	78*300 mm	20 ng/mL
SS181-3-U	THC Rapid Test Uncut Sheet	Urine	60*300 mm	40 ng/mL
SS181-4-U	THC Rapid Test Uncut Sheet	Urine	78*300 mm	40 ng/mL
SS181-5-U	THC Rapid Test Uncut Sheet	Urine	60*300 mm	50 ng/mL
SS181-6-U	THC Rapid Test Uncut Sheet	Urine	78*300 mm	50 ng/mL
SS181-1-S	THC Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS186-1-U	MDMA Rapid Test Uncut Sheet	Urine	60*300 mm	500 ng/mL
SS186-2-U	MDMA Rapid Test Uncut Sheet	Urine	78*300 mm	500 ng/mL
SS186-1-S	MDMA Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS195-1-U	TRA Rapid Test Uncut Sheet	Urine	60*300 mm	100 ng/mL
SS195-2-U	TRA Rapid Test Uncut Sheet	Urine	78*300 mm	100 ng/mL
SS191-1-U	FTY Rapid Test Uncut Sheet	Urine	60*300 mm	25 ng/mL
SS191-2-U	FTY Rapid Test Uncut Sheet	Urine	78*300 mm	25 ng/mL
SS191-1-S	FTY Rapid Test Uncut Sheet	Saliva	60*300 mm	10 ng/mL
SS191-3-U	NFTY Rapid Test Uncut Sheet	Urine	60*300 mm	10 ng/mL
SS191-4-U	NFTY Rapid Test Uncut Sheet	Urine	78*300 mm	10 ng/mL
SS191-3-S	NFTY Rapid Test Uncut Sheet	Saliva	60*300 mm	10 ng/mL
SS182-1-U	MTD Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS182-2-U	MTD Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS182-1-S	MTD Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS193-1-U	EDDP Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS193-2-U	EDDP Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS194-1-U	K2 Rapid Test Uncut Sheet	Urine	60*300 mm	50 ng/mL
SS194-2-U	K2 Rapid Test Uncut Sheet	Urine	78*300 mm	50 ng/mL
SS194-1-S	K2 Rapid Test Uncut Sheet	Saliva	60*300 mm	25 ng/mL

Drug Abuse

Cat. No.	Uncut Sheet	Specimen	Specification	Cut Off
SS189-3-U	OPI Rapid Test Uncut Sheet	Urine	60*300 mm	2000 ng/mL
SS189-4-U	OPI Rapid Test Uncut Sheet	Urine	78*300 mm	2000 ng/mL
SS189-2-S	OPI Rapid Test Uncut Sheet	Saliva	60*300 mm	40 ng/mL
SS179-1-U	COC Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS179-2-U	COC Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS179-1-S	COC Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS180-1-U	AMP Rapid Test Uncut Sheet	Urine	60*300 mm	1000 ng/mL
SS180-2-U	AMP Rapid Test Uncut Sheet	Urine	78*300 mm	1000 ng/mL
SS180-1-S	AMP Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS178-1-U	BZO Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS178-2-U	BZO Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS178-1-S	BZO Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS190-1-U	BAR Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS190-2-U	BAR Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS190-1-S	BAR Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS183-1-U	TCA Rapid Test Uncut Sheet	Urine	60*300 mm	1000 ng/mL
SS183-2-U	TCA Rapid Test Uncut Sheet	Urine	78*300 mm	1000 ng/mL
SS184-1-U	BUP Rapid Test Uncut Sheet	Urine	60*300 mm	10 ng/mL
SS184-2-U	BUP Rapid Test Uncut Sheet	Urine	78*300 mm	10 ng/mL
SS184-1-S	BUP Rapid Test Uncut Sheet	Saliva	60*300 mm	10 ng/mL
SS245-1-U	COT Rapid Test Uncut Sheet	Urine	60*300 mm	100 ng/mL
SS245-2-U	COT Rapid Test Uncut Sheet	Urine	78*300 mm	100 ng/mL
SS154-1-U	ETG Rapid Test Uncut Sheet	Urine	60*300 mm	1000 ng/mL
SS154-2-U	ETG Rapid Test Uncut Sheet	Urine	78*300 mm	1000 ng/mL
SS246-1-U	MQL Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS246-2-U	MQL Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS188-1-U	OXY Rapid Test Uncut Sheet	Urine	60*300 mm	100 ng/mL
SS188-2-U	OXY Rapid Test Uncut Sheet	Urine	78*300 mm	100 ng/mL
SS188-1-S	OXY Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS185-1-U	PCP Rapid Test Uncut Sheet	Urine	60*300 mm	25 ng/mL
SS185-2-U	PCP Rapid Test Uncut Sheet	Urine	78*300 mm	25 ng/mL
SS185-1-S	PCP Rapid Test Uncut Sheet	Saliva	60*300 mm	25 ng/mL
SS187-1-U	PPX Rapid Test Uncut Sheet	Urine	60*300 mm	300 ng/mL
SS187-2-U	PPX Rapid Test Uncut Sheet	Urine	78*300 mm	300 ng/mL
SS187-1-S	PPX Rapid Test Uncut Sheet	Saliva	60*300 mm	50 ng/mL
SS189-5-U	6-AM Rapid Test Uncut Sheet	Urine	60*300 mm	10 ng/mL
SS189-6-U	6-AM Rapid Test Uncut Sheet	Urine	78*300 mm	10 ng/mL
SS189-3-S	6-AM Rapid Test Uncut Sheet	Saliva	60*300 mm	10 ng/mL

Order Information:

HomePage: <http://heavybio.com/>

E-mail: sales@heavybio.com

Address: Building 47, Baotian Industrial Zone, Chentian Community, Xixiang Street, Bao'an District, Shenzhen City, Guangdong Province, China